BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014921	Date: 08-JAN-2020	Proofread No.: 1
Designer: CDS	Page: 1 of 16	
Colour Reference:	BLACK	

prismaflex

M60 SET / M100 SET / M150 SET

AN69 HF
MEMBRANE

Instructions for use
使用説明

REF 115305
REF 115306
REF 115307

Made in France

Baxter

CE 2797

Date of Revision: 2019/12/01

Blood warmer connection (blue)
血液加溫器連接 (藍色)

BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014921	Date: 08-JAN-2020	Proofread No.: 1
Designer: CDS	Page: 3 of 16	
Colour Reference:	BLACK	

DEFINITION OF SYMBOLS

		ENGLISH		繁體中文	
SYMBOL GRAPHIC & REF NUM / 符號圖形和 參考編號	TITLE AND NUMBER OF STANDARD / 標準名 稱和編號	SYMBOL TITLE	SYMBOL DESCRIPTION (EXPLANATORY TEXT)	符號名稱	符號描述（解釋性文字）
5.1.1 	ISO 15223-1 ¹	Manufacturer	Indicates the medical device manufacturer, i.e. the natural or legal person with responsibility for the design, manufacture, packaging and labelling of a device before it is placed on the market under his own name.	製造廠	表示醫療器材製造廠，即在其名義將器材 上市前負責器材設計、生產、包裝和標籤的自然 人或法人。
5.1.3 	ISO 15223-1 ¹	Date of manufacture	Indicates the date when the medical device was manufactured. Format should be YYYY-MM-DD.	製造日期	表示生產該醫療器材的日 期。格式應為年-月-日。
5.1.4 	ISO 15223-1 ¹	Use-by date	Indicates the date after which the medical device is not to be used. Format should be YYYY-MM-DD	保存期限	表示在該日期後不應使用 該醫療器材。格式應為 年-月-日
5.1.5 	ISO 15223-1 ¹	Batch code	Indicates the manufacturer's batch code. Synonyms are "lot number" and "batch number"	批次代碼	表示製造廠的批次代碼。同義字為“批號”和“批次號”
5.1.6 	ISO 15223-1 ¹	Catalogue number	Indicates the manufacturer's catalog number. Synonyms are "reference number" and "reorder number"	目錄編號	表示製造廠的目錄編號。同義字為“參考編號”和“再訂購編號”
5.2.3 	ISO 15223-1 ¹	Sterile fluid path	Indicates the presence of a sterile fluid path within the medical device in cases where other parts of the medical device, including the exterior, might not be supplied sterile. The method of sterilization shall be indicated in the empty box, as appropriate.	無菌液體通路	表示當該醫療器材其它部 分（包括外部）不能以無 菌形式供應情況下，該醫 療器材的內部為無菌液體 通路。應將滅菌方法標示 在外盒。
5.2.8 	ISO 15223-1 ¹	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.	如果包裝破損，請勿使用	表示如果包裝已破損或打 開，不應使用醫療器材
5.3.1 	ISO 15223-1 ¹	Fragile, handle with care	Indicates a medical device that can be broken or damaged if not handled carefully.	易碎、小心輕放	表示如果不注意操作， 可能會打破或損壞醫療 器材。
5.3.4 	ISO 15223-1 ¹	Keep dry OR Keep away from rain	Indicates a medical device that needs to be protected from moisture. OR Indicates a medical device that needs to be kept away from rain.	保持乾燥 或 遠離雨水	表示醫療器材需要防潮。 或 表示醫療器材需要遠離 雨水。
5.3.7 	ISO 15223-1 ¹	Temperature Limit	Indicates the temperature limits to which the medical device can be safely exposed.	溫度限值	表示該醫療器材可安全地 暴露於該溫度限值內。
5.4.2 	ISO 15223-1 ¹	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure. Synonyms for "Do no re-use" are "single use" and "use only once".	請勿重複使用	表示醫療器材僅供單次使 用，或用於單一病人的單 次治療。“ 請勿重複使 用”的同義字為“單次使 用”和“僅使用一次”。
5.4.4 	ISO 15223-1 ¹	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.	注意	表示使用者需要參閱使用 說明、了解重要警訊，如 因為各種原 因不能標示在 醫療器材機體上的各種警 告和預防措施。
	NA, European countries only / NA, pays européens seulement / NA, sólo países europeos	Green Dot Symbol	To indicate the manufacturer of the product contributes to the cost of recovery and recycling	綠點符號	表示該產品製造廠落實回 收和重複利用
0623 	ISO 7000 ²	This way up	Indicates the correct, upright position of the package	此面朝上	表示包裝的正確、朝上位置
5.4.3 	ISO 15223-1 ¹	Consult instructions for use	Indicates the need for the user to consult the instructions for use. Synonym for "Consult instructions for use" is "Consult operating instructions".	參閱使用說明	表示使用者需要參閱使用 說明。“參閱使用說明” 的同義字為 “參閱操作 手冊”。

¹ISO 15223-1:2016 Medical Devices - Symbols to be used with medical device labels, labelling and information to be supplied Part 1: General Requirements / 醫療器材 - 要提供與醫療器材標籤配合使用的符號、標記和資訊第 1 部分：基本要求
²ISO 7000 Graphical symbols for use on equipment - Registered symbols / 供設備上使用的圖形符號 - 已註冊符號

BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014921	Date: 08-JAN-2020	Proofread No.: 1
Designer: CDS	Page: 4 of 16	
Colour Reference:	BLACK	

eng	English	5
chi	繁體中文	9
	Installation guide	15

BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014921	Date: 08-JAN-2020	Proofread No.: 1
Designer: CDS	Page: 5 of 16	
Colour Reference:	BLACK	

ENGLISH

The **PRISMAFLEX M60/M100/M150 Set** is manufactured by **GAMBRO Industries, 7 avenue Lionel Terray, BP 126, 69883 MEYZIEU CEDEX, FRANCE.**

⚠ Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

DEFINITION OF EXPRESSIONS USED
IN THIS MANUAL

In this document :

⚠ **“Warning”** is used to alert the user/operator not to take a certain action which, if taken, can cause a potential hazard and result in a serious adverse reaction, injury or death.

⚠ **“Caution”** is used to alert the user/operator to take a certain action to protect against a potential hazard, which, if ignored, could have an adverse effect on the patient or on the device.

“Note” is used as a reminder to the user/operator on normal treatment activity and on what is a suitable action in a particular situation.

SCUF: Slow Continuous UltraFiltration.
CVVH: Continuous Veno-Venous Hemofiltration.
CVVHD: Continuous Veno-Venous HemoDialysis.
CVVHDF: Continuous Veno-Venous HemoDiaFiltration.
Predilution : addition of replacement fluid to the blood stream **upstream** to the filter.
Postdilution: addition of replacement fluid to the blood stream **downstream** to the filter.

PRODUCT DESCRIPTION

- The **PRISMAFLEX M60/M100/M150 Set** is a disposable, extracorporeal circuit for use with the **PRISMAFLEX** system.
- The **PRISMAFLEX M60/M100/M150 Set** consists of a AN69 HF hollow fiber hemofilter/dialyzer* and tubing lines.
- This filter is permanently connected to a blood access line (red-striped), a blood return line (blue-striped), a dialysate inlet line (green-striped) and an effluent outlet line (yellow-striped).
- The other lines of the set include:
 - a replacement solution line (purple-striped),
 - a pre blood pump line (white striped),
 - an anticoagulant line (syringe).
- The configuration of the **PRISMAFLEX** set allows the following uses depending on the configuration of the automated clamps on the machine:
 - purple circuit: replacement in pre or post dilution (CVVH and CVVHDF),
 - green circuit:
 - * dialysate in CVVHD and CVVHDF,
 - * replacement in post dilution in CVVH.
- The pre blood pump line allows the addition of infusion solution close to the end of the patient access line and before the blood pump. This can be used as an additional pre-dilution infusion to the replacement circuit.
- The **PRISMAFLEX** Set is provided with a specific small volume deaeration chamber in which blood does not appear to mix with the replacement liquid the majority of the time; this is a normal operation of the device.
- A 5-liter bag is provided to be connected to the end of the blood return line to initially collect priming solution, during priming. Then, during treatment, this bag is used to collect ultrafiltrate and/or used dialysate (connection at effluent line). Other sterile 5 and 9 liter bags and sterile, non pyrogenic spikes can be ordered separately.
- All line connectors are compatible with the ISO 594/1 & 2 international standards concerning conical fittings.
- The fluid pathways of the **PRISMAFLEX** Set are guaranteed sterile and non pyrogenic.
- The **PRISMAFLEX** Set is sterilized by ethylene oxide (EtO). Deaeration is such that EtO residuals comply with the ones described in ISO 10993-7.
- Expiration date: please refer to product label.

* In this document the hemofilter/dialyzer will be referred to as “filter”.

INTENDED USE / INDICATIONS

The **PRISMAFLEX** Set is indicated for use only with the **PRISMAFLEX** control unit in providing continuous fluid management and renal replacement therapies. The system is intended for patients who have acute renal failure, fluid overload, or both.

This set is intended for use in the following veno-venous therapies: SCUF; CVVH; CVVHD; CVVHDF.

All treatments administered via the **PRISMAFLEX** Set must be prescribed by a physician. The size, weight, state of uremia, cardiac status, and general physical condition of the patient must be carefully evaluated by the prescribing physician before each treatment.

CONTRAINDICATIONS

There are no known contraindications to continuous renal replacement therapies.

CAUTIONS AND WARNINGS

Note: additional warnings and cautions pertaining to the **PRISMAFLEX** system are included in the **PRISMAFLEX** control unit operator’s manual.

⚠ Cautions

It is recommended that particular attention must be paid with respect to extra-corporal blood volume (see General Characteristics).

The pediatric use of the **PRISMAFLEX M60 Set** should be restricted to children with a body weight greater than 11 kg (24lb) with respect to extracorporeal blood volume. Refer to **PRISMAFLEX** operator’s manual for any additional weight restrictions that might apply.

The **PRISMAFLEX M100 Set** and **M150 Sets** should be restricted to patients with a body weight greater than 30kg (66lb).

1. Carefully read these instructions for use and the **PRISMAFLEX** control unit operator’s manual before using this product.
2. Store the **PRISMAFLEX** Set in a dry place, between 0° C (32° F) and 30° C (86° F).
3. Some solvents and other chemicals, if used in contact with the filter, could damage the set. No chemical of this type should be used without permission of the manufacturer. The following are especially forbidden:
 - a) halogenated aromatic and aliphatic solvents,
 - b) ketonic solvents.
4. To prevent contamination, this **PRISMAFLEX** Set must be used as soon as its packaging and sterilization caps are removed.
5. Do not use this set if the packaging is damaged, if the sterilization caps are missing or loose, or if any of the lines in the set are kinked.
6. Do not try to remove the filter from the cartridge plate.
7. Destroy this set after single use, using aseptic technique for potentially contaminated equipment. Do not re-sterilize. The **PRISMAFLEX M** is intended for single use only. Re-using the **PRISMAFLEX M** may cause serious damage to the product resulting in patient injury or death.
8. Use aseptic techniques when handling all blood and fluid lines in the set.
- 9a. Use only prescribed dialysate and replacement solutions with the **PRISMAFLEX** system. These solutions must have a density similar to that of saline solutions (close to 1) in order to avoid errors in the volumes used for fluid exchange.
- 9b. In CVVHD and CVVHDF modes, it is recommended to use only sterile bagged dialysate.
- 9c. In CVVH modes and CVVHDF, if a commercially available replacement solution is used, it must be labeled as intended for intravenous injection.
10. Connect the **PRISMAFLEX** Set to a patient via venous blood access and return devices. A double-lumen venous catheter is the recommended blood access device; however, two single-lumen venous catheters can also be used. There are 3 possible accesses for **PRISMAFLEX** system therapies: subclavian, jugular or femoral vein.
11. During priming and operation, observe closely for leakage at joints and connections within the set, notably the bags. Leakage can cause blood loss or air embolism. If leakage cannot be stopped by tightening the connections, replace the set.

BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014921	Date: 08-JAN-2020	Proofread No.: 1
Designer: CDS	Page: 6 of 16	
Colour Reference:	BLACK	

12.

Before connecting the blood return line to the patient, check for absence of air between the segment of line inserted in the air detector and the patient-end of the return line.

If air is present in this part of the return line, connect the access line to the patient and start the blood pump while leaving the return line connected to the collection bag. Purge the air present in the end-part of the return line, then stop the blood pump. Disconnect the return line from the collection bag, and connect it to the patient.

If the amount of air in the blood circuit is too large, reprime the circuit completely before patient connection.
13.

After priming is complete, do not remove the pressure pods from the pressure sensor housings. If pods are removed, the set must be changed or the Diaphragm Reposition procedure performed (refer to **PRISMAFLEX** control unit operator's manual).
14.

If the patient is not immediately connected to the **PRISMAFLEX** Set after priming is complete, flush the set with at least 1000 mL priming solution [saline or alkaline solution (pH ≥ 7.3) with heparin added] prior to connecting the patient. This requires use of a new bag of priming solution.
15.

Use a 21-gauge or smaller needle to obtain blood/fluid samples or remove trapped air from the **PRISMAFLEX** Set. Use of larger needles can cause holes in the sample sites, resulting in blood loss or air embolism.
16.

The **PRISMAFLEX** control unit may not be able to detect disconnections of the set from the patient's catheter. Carefully observe the set and all operations while using the **PRISMAFLEX** system for a patient treatment.
17.

Due to the nature of use of the **PRISMAFLEX** Set (low blood flow rates, extended treatment time, and other special factors), the possibility for coagulation within the blood flowpath is substantially enhanced. Give careful attention to the possible medical hazards associated with coagulation of the blood flowpath and comply with the minimum blood flow rates specifications of each filter (see the "Filter Operating Specifications" section).
18.

In the case of patients who pose a high risk of hemorrhaging it is recommended not to add heparin to the priming solution.
19.

Filter performance specifications require a minimum blood flow rate, specific to each filter, to avoid risk of hemoconcentration (see "Filter Operating Specifications" section).
20.

During use, closely monitor the patient's clotting parameters, especially when increasing the amount of anticoagulant delivered or after changing the anticoagulant syringe.
21.

To assure proper anticoagulant flow control, only use the syringes listed in the operator's manual. The use of non-recommended syringes can be a hazard for the patient. Particularly if there is no Luer-lock on the syringe, the seal between the syringe and the heparin line can no longer be guaranteed.
22.

When not using the pre blood pump infusion circuit, it is recommended to clamp this circuit close to its connection to the access line; this will prevent the sedimentation of blood into the pre blood infusion line.
23.

Always inspect the blood flowpath for signs of clotting before returning the blood in the set to the patient. If clotting is suspected, do not return the blood to the patient.
24.

In case of re-circulation mode, the set must be replaced if the maximum re-circulation time is exceeded; refer to the **PRISMAFLEX** operator's manual for more information.

In case of poor blood return, the set must be replaced.

In all cases, it is essential to reprime the set with fresh saline immediately before patient connection.
25.

The **PRISMAFLEX** set offers a specific design of the deaeration chamber which aims to trap air before blood is returned to the patient.
26.

The **PRISMAFLEX** Set is not designed for a heater to be connected to the replacement solution line. A heater generates air bubbles which collect in the return deaeration chamber. Therefore, it is recommended not to use a heater on the replacement solution line.
27.

The **PRISMAFLEX** set is not designed for a heater to be connected to the dialysate solution line. A heater generates air bubbles which collect in the filtrate/dialysate compartment of the filter and decrease diffusive performance of the device. Therefore it is recommended not to use a heater on the dialysate solution line.

In general terms, introduction of air in the dialysate circuit should be minimised throughout the course of the treatment, especially when replacing dialysate solution bags.
28.

Hypothermia must be monitored in all CRRT treatments, and special attention should be paid when increasing exchange volumes above 2 L/h ; it may be necessary to warm the patient because of hypothermia.
29.

The deaeration chamber line (blue-striped) is equipped with a Luer-lock connection near the deaeration chamber.

This connector is intended to join the extension line of a blood warmer. Refer to the specific Instructions for Use and strictly follow detailed instructions for set up of this line. Do not use this connection for any other purpose.
30.

Do not attach/connect the extension line of a blood heater to the return line downstream of the air detector. The **PRISMAFLEX** system cannot detect air introduced in the line downstream of the air detector.

⚠

Warnings

1.

The use of operating procedures other than those published by the manufacturer or the use of accessory devices not recommended by the manufacturer can result in patient injury or death.
2.

Use only **PRISMAFLEX** Sets with the **PRISMAFLEX** control unit. The use of non-**PRISMAFLEX** Sets can result in patient injury or death.
3.

Should acute allergic reactions (first-use syndrome) occur in patients receiving treatment, **immediately stop the treatment** and administer appropriate intervention. Pay special attention to patients receiving ACE inhibitors and/or having already shown similar allergic reactions (see "Hypersensitivity Reactions" section).
4.

Do not allow air to enter the blood compartment of the filter after priming is started. If a large amount of air enters, the set must be replaced.
5.

Since drugs can pass through the membrane of the filter, the dosage of associated drug treatments must be adjusted for patients on continuous renal replacement therapy.
6.

To assure adequate filter performance, it is recommended that the set be changed every 24 hours of use. However, the set must be changed after 3 days (72 hours) and/or the maximum process volume of blood (780 L) whichever occurs first. Continued use beyond these limits (either 72 hours or 780 L) could result in rupture of the pump segments, with risk of patient injury or death.
7.

Use only drugs compatible with plastics listed in the specifications section. Some plastics can be incompatible with drugs when in contact with solutions with pH > 10.

SPECIFICATIONS

See Tables at end of document.

SET MATERIALS

- AN69 HF hollow fiber

: Acrylonitrile and sodium methallyl sulfonate copolymer
- Housing and headers

: Polycarbonate
- Potting compound

: Polyurethane
- Tubing material

: Plasticized polyvinyl chloride (PVC)
- Cartridge

: PETG

Note: the following information is available from the manufacturer upon request:

- information about test methods used to obtain performance characteristics,
- the number and range of particles in the effluent from the dialyzer prepared as recommended for clinical use,
- the types and amounts of residue from the sterilization process.

The **PRISMAFLEX** Set is not made with rubber natural latex.

INSTRUCTIONS FOR USE

Note: use the set by following the detailed on-line instructions provided by the **PRISMAFLEX** control unit. Additional information is available in the **PRISMAFLEX** control unit operator's manual.

Note: a TMP > 40 kPa (300 mmHg) does not allow a higher ultrafiltration.

Perform the following procedures when the appropriate instructions appear on the display of the **PRISMAFLEX** control unit.

Load Set

1.

Remove the set from the packaging support. Holding the filter vertically (so that the label is the right way up), carefully snap the set cartridge into the cartridge carrier (center of front panel).
2.

Attach the 3 pressure pods to their proper pressure housings. Press effluent line into blood leak detector; snap discharger ring into its guide.
3.

Temporarily hang access/effluent Y line on priming hook.
4.

Place deaeration chamber in its holder; attach chamber monitor line to return pressure port.
5.

Insert return line into air detector and return line clamp.
6.

Connect return line to effluent bag.
7.

Open effluent scale; hang collection/effluent bag. Close scale.

Prepare and Connect Solutions

1.

Hang bag of priming solution [saline or alkaline solution (pH ≥ 7.3) with added 5000IU heparin/liter] on priming hook (left corner hook top of front panel). Connect access (red)/effluent (yellow) Y-line to priming solution bag.
2.

If required, connect PBP line (white) to pre blood pump (PBP) bag; hang bag on its scale.
3.

Hang replacement solution (CVVH, CVVHDF) on purple scale hook. Connect replacement solution line (purple) to replacement solution bag.

BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014921	Date: 08-JAN-2020	Proofread No.: 1
Designer: CDS	Page: 7 of 16	
Colour Reference:	BLACK	

4. In CVVHD/CVVHDF hang dialysate on green scale hook. Connect dialysate line (green). In CVVH hang replacement solution on green scale hook (postdilution replacement). Connect green striped line to bag.
Note: see Caution no. 10 a, b, c.
5. Connect anticoagulant line to filled anticoagulant syringe. Install syringe in pump (see Help).
6. Unclamp any clamped lines. **Verify all connections are secure.** Press PRIME to start automatic priming.

Prime Set

Note: see Cautions no. 12 through 15, cautions no. 27, 28 and 31, and Warning no. 4.

Priming includes multiple self-tests and takes approximately 10 minutes. After the cycle is complete:

1. Examine set carefully to be certain all connections are secure, all lines are unobstructed, and there are no leaks in the tubing.
2. Leave priming solution and prime collection bags attached until ready to connect patient.
3. Continue chosen treatment by following the instructions on the display of the **PRISMAFLEX** control unit.

The **PRISMAFLEX** Set must be carefully deaerated.

Anticoagulation Considerations

Note: see Cautions no. 18 through 22 and Warning no. 6.

Initiate anticoagulation of the blood flowpath, as prescribed by the physician. During use, monitor the patient’s clotting parameters; adjust the anticoagulation settings on the **PRISMAFLEX** control unit, according to the physician’s prescription. If prescribed, do not forget to infuse a loading dose of anticoagulant immediately after patient connection.
Anticoagulation plays an important part in extending filter life by retarding plugging and clotting.

Change Set Procedure

To remove this set, load a new set and continue with present treatment: Press “STOP” from the Status screen, then press “CHANGE SET” and follow the on-line instructions.
Note: operator can return blood to the patient prior to disconnecting, if desired (see Caution no. 24).

Re-circulation procedure

To set-up the re-circulation of the circuit:
Press “STOP” on the status screen; then press “RECIRC” and follow the on-line instructions.

Note: the operator must return blood present in the set to the patient, then disconnect the patient and circulate a sterile saline solution in the blood circuit of the set. Once the treatment can be re-started, the set must be re-primed and re-rinsed with a sterile saline solution before reconnecting the patient (see cautions n° 24 and 25).

End Treatment Procedure

To end the present treatment and remove this set :
Press “STOP” from the Status screen, then press “END TREATMENT” and follow the on-line instructions.
Note: operator can return blood to the patient prior to disconnecting, if desired (see Caution no. 24).

MANUAL TERMINATION

Manual termination may be necessary due to power loss or an alarm of the **PRISMAFLEX** control unit. The alarm screen tells the operator if a manual termination is required.
Note: the following instructions are also found in “Troubleshooting” in the **PRISMAFLEX** control unit operator’s manual.

A. With Blood Return

Note: see Caution no. 24.

Note: a sterile spike connector may be required.

1. Turn off the power. Clamp the access line (red-striped) and disconnect from the patient. Attach the access line to a 1-liter bag of sterile saline (use spike connector if needed). Unclamp the access line.
2. Press the return clamp button (left side of the return line clamp assembly) and hold in the “in” position. With the other hand, remove the return line (blue-striped) from the return line clamp.
3. Visually check the fluid level in the deaeration chamber. If the level is insufficient:
- disconnect the deaeration chamber service line from the return pressure port on the **PRISMAFLEX** machine (the level will automatically rise in the deaeration chamber),
- reconnect the line once the correct fluid level is reached.
4. Remove the pump crank from its holder on the rear panel. Insert crank into the rotor of the blood pump and turn clockwise until sufficient blood is returned to the patient.
 **Warning:** the alarm system is disabled. Visually check for air in the blood return line until the patient is disconnected.
5. Clamp the return line (blue-striped) and disconnect from the patient. Clamp lines to all bags.
6. Press the two clips of the cartridge carrier to release the cartridge. Starting with the peristaltic pump, insert the pump crank into the rotor and turn each pump counterclockwise.
7. When the pump segments are free, grasp the cartridge and pull out to disengage the lines from the pinch valves. Take the set off the control unit and discard as usual.

Note: remaining solutions may be used with a new set, if desired.

 **Warning:** ensure patient is disconnected from set before removing set from control unit.

B. Without Blood Return

Note: the patient will lose the blood contained in the blood flowpath during a manual termination without blood return.

1. Turn off the power. Clamp the access line (red-striped) and return line (blue-striped) and disconnect from the patient.
2. Clamp lines to all bags.
3. Press the two clips of the cartridge carrier to release the cartridge. Starting with the blood pump, insert the pump crank into the rotor and turn each pump counterclockwise.
4. When the pump segments are free, grasp the cartridge and pull out to disengage the lines from the pinch valves. Take the set off the control unit and discard as usual.

SPECIAL PROCEDURES IN CASE OF COMPLICATION

Filter Membrane Blood Leaks

Blood leaks through the filter membrane are automatically detected by the **PRISMAFLEX** control unit alarm system. A warning alarm is generated and blood loss is limited by immediate stoppage of all pumps.
To return blood to the patient, press STOP from the alarm screen, then press CHANGE SET from the Stop screen and follow the on-line instructions.

External Blood Leaks

Note: see Cautions no. 16, 17 and 22.

External blood leakage may not be immediately identified by monitoring equipment and could result in significant blood loss. Check the filter and all connections of the disposable tubings during treatment to minimize the risk of leakage. If an external blood leakage is observed, immediately stop the blood pump. Initiate corrective action by securing connections or replacing the **PRISMAFLEX** Set.
If necessary, administer adequate replacement solution to the patient to compensate for blood loss.

Hypersensitivity Reactions

Note: see Warning no. 3.

Should acute allergic reactions (first use syndrome) occur within the first few minutes of the treatment, it is important to react immediately by discontinuing the session and administering appropriate treatment.

BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014921	Date: 08-JAN-2020	Proofread No.: 1
Designer: CDS	Page: 8 of 16	
Colour Reference:	BLACK	

Patients receiving angiotensin converting enzyme (ACE) inhibitors as medication can develop, within the first few minutes of a treatment, symptoms similar to acute allergic reactions i.e bronchospasm, edema of airways or larynx, dyspnea, angioedema, urticaria, nausea, vomiting, diarrhea, respiratory arrest, abdominal cramping, hypotension, hypovolemic shock and death.

However, for these patients, administration of antihistamines often does not alleviate the symptoms. In this case, treatment must be stopped and a more aggressive first-line therapy for an anaphylactoid reaction should be initiated immediately after the onset of symptoms.

Therefore, special attention ought to be paid to patients receiving ACE inhibitors and/or having already shown similar reactions.

WARRANTY AND LIMITATION OF LIABILITY

a) The manufacturer warrants that the **PRISMAFLEX** Set has been manufactured in accordance with its specifications and in compliance with good manufacturing practices, other applicable industry standards and regulatory requirements.

If provided with the lot/serial number of the defective product, the manufacturer will, by replacement or credit, remedy manufacturing defects in the **PRISMAFLEX** Set becoming apparent before the expiration date.

b) The warranty under paragraph a) above is in lieu of, and to the exclusion of, any other warranty, whether written or oral, express or implied, statutory or otherwise, and there are no warranties of merchantability or other warranties, which extend beyond those described in paragraph a) above. The remedy set out above for manufacturing defects is the sole remedy available to any person due to defects in the **PRISMAFLEX** Set and the manufacturer shall not be liable for any consequential or incidental loss, damage, injury or expense arising directly or indirectly from the use of the **PRISMAFLEX** Set, whether as a result of any defect therein or otherwise.

c) The manufacturer shall not be liable for any misuse, improper handling, non-compliance with warnings and instructions, damage arising from events after the manufacturer's release of the **PRISMAFLEX** Set, failure or omission to inspect the **PRISMAFLEX** Set before use in order to ensure that the **PRISMAFLEX** Set is in proper condition, or any warranty given by independent distributors or dealers.

d) The manufacturer is GAMBRO Industries, 7 avenue Lionel Terray, BP 126, 69883 MEYZIEU CEDEX, FRANCE.

INDIA

- **Name and address of importer:** Baxter India Pvt. Ltd., Gala No. 1 to 6, Building No. C15, Ground Floor, Shree Arihant complex, Kalher, Bhiwandi, Dist. Thane, Tal: Bhiwandi (Thane-Zone 5) 421302
- Name, address, telephone no., email address (in case of consumer complaint):

Baxter India Pvt. Ltd.

5th Floor, Block A, Building 9,

DLF Phase III, Cyber City, Gurgaon- 122002.

Consumer Care No.: 0124-4603200

Consumer Care email ID: Customerservice_SHS_India@baxter.com
- Import License No.: MD-1591-2310

BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014921	Date: 08-JAN-2020	Proofread No.: 1
Designer: CDS	Page: 9 of 16	
Colour Reference:	BLACK	

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BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014921		Date: 08-JAN-2020 Proofread No.: 1
Designer: CDS		Page: 12 of 16
Colour Reference: <div>BLACK</div>		

FILTER DATA	過濾器資料	PRISMAFLEX		
		M60 SET	M100 SET	M150 SET
NOMINAL PHYSICAL CHARACTERISTICS Effective surface area Fiber internal diameter (wet) Fiber wall thickness	標稱物理特性 有效表面面積 纖維內徑 (潮濕) 纖維壁厚	0.60 m²	0.90 m²	1.50 m²
		240 µm		
		50 µm		
IN VITRO PERFORMANCES ≠ Blood priming volume	體外效能 ≠ 血液預充量	42 mL ± 10%	66 mL ± 10%	105 mL ± 10%
Blood pressure drop (post dilution) (bovine blood, Htc****32%, Cp***** 60 g/L, 37°C) QB** = 100 mL/min, QUF*** = 1 L/h QB** = 180 mL/min, QUF*** = 2 L/h QB** = 250 mL/min, QUF*** = 2 L/h QB** = 400 mL/min, QUF*** = 2 L/h QB** = 450 mL/min, QUF*** = 2L/h	血壓壓力降 (後稀釋) (牛血, Htc****32%, Cp***** 60 g/L, 37°C) QB** = 100 mL/min, QUF*** = 1 L/h QB** = 200 mL/min, QUF*** = 2 L/h QB** = 300 mL/min, QUF*** = 2 L/h QB** = 400 mL/min, QUF*** = 2 L/h QB** = 450 mL/min, QUF*** = 2 L/h	± 20 % 47 mmHg 91 mmHg	± 20 % 31 mmHg 60 mmHg 73 mmHg 105 mmHg	± 20 % 20 mmHg -- 51 mmHg 64 mmHg 70 mmHg
Sieving coefficient (bovine plasma, Cp 60 g/l, 37°C) QB** = 100 mL/min, QUF*** = 20 mL/min - Urea - Creatinine - Vitamin B12 - Inulin - Myoglobin - Albumin	篩濾係數 (牛血漿, Cp 60 g/L, 溫度 37° C) QB** = 100 mL/min, QUF*** = 20 mL/min - 尿素 - 肌酸酐 - 維生素 B12 - 菊糖 - 肌紅蛋白 - 白蛋白	1 1 1 0.95 0.55 <0.01		

≠ Typical mean values obtained from laboratory testing of post-sterilization sample lots. Results may vary depending on patient and clinical conditions.
由實驗室測試滅菌後樣品所得之平均值。視病人及臨床狀況而定，結果可能不同。

GENERAL DATA / 一般資料

	PRISMAFLEX		
	M60 SET	M100 SET	M150 SET
Weight / 重量	780 g	800 g	860 g
Overall dimensions / 整體尺寸			
lenght / 長度	27 cm	27 cm	27 cm
width / 寬度	22 cm	22 cm	22 cm
height / 高度	9 cm	9 cm	9 cm
Blood volume in set / 套組中的血液容積	93 mL	152 mL	189 mL




FILTER OPERATING SPECIFICATIONS

過濾器操作規格	PRISMAFLEX		
	M60 SET	M100 SET	M150 SET
Maximum TMP* / 最高	450 mmHg 60 kPa		
Maximum blood pressure / 最高血壓	500 mmHg 66,6 kPa		
Minimum blood flow rate / 最低血液流速	50 mL/min	75 mL/min	100 mL/min

		PRISMAFLEX M60 SET				PRISMAFLEX M100 SET				PRISMAFLEX M150 SET			
		QB/QS** = 100 mL/min QUF*** = 0 mL/min				QB/QS**= 150 mL/min QUF*** = 0 mL/min				QB/QS**= 200 mL/min QUF*** = 0 mL/min			
QD*****		(L/h)	1	2,5	4	1	2,5	4	8	1	2,5	4	8
		(mL/min)	17	42	67	17	42	67	133	17	42	67	133
Urea / 尿素		(±10%)	17	39	54	17	41	63	95	17	42	66	117
Vit B12 / 維生素 B12		(±20%)	14	23	28	16	30	37	45	17	37	49	64
Inulin / 菊糖		(±20%)	12	17	19	14	23	26	30	16	31	37	45

CVVHD CLEARANCES (Continuous veno-venous hemodialysis) Clearances versus inlet dialysate flow rate (37°C)
CVVHD 廓清率 (連續性靜脈-靜脈血液透析) 廓清率比入口透析液流速 (37°C)

- * Transmembrane pressure / 跨膜壓
- ** Arterial blood flow rate / 動脈血液流速
- *** Ultrafiltration flow rate (on PRISMA system, the ultrafiltration flow rate = fluid removal flow rate + replacement flow rate) / 超過濾流速 (在 PRISMA 系統上，超過濾流速 = 液體清除流速 + 置換液流速) /
- **** Hematocrit / 血容積比
- ***** Protein concentration / 蛋白濃度
- ***** Dialysate flow rate / 透析液流率

“In vitro” ultrafiltration with blood (values ± 15%).
(Bovine blood at 37°C, Hct 32 %, Protein concentration 60 g/L).
Ultrafiltration is controlled by the PRISMAFLEX System and is independent of the ultrafiltration coefficient (KUF).
「體外」血液超過濾 (數值 ±15 %).
(37°C 牛血，Hct 32 %，蛋白質濃度 60 g/L).
超過濾係由 PRISMAFLEX 系統控制，與超過濾係數 (KUF) 無關

BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014921	Date: 08-JAN-2020	Proofread No.: 1
Designer: CDS	Page: 14 of 16	
Colour Reference:	BLACK	

Prismatherm




SP420 line connexion




INSTALL SET DIRECTLY ONTO

prismaflex

MONITOR

Patented ready to use PrismaFlex set packaging.